


Product Preparation and Recommendations



With the same strength as the ZirFit® Solid Zirconia, ZirFit® HT offers high translucency with very high esthetic. ZirFit® solutions are all ceramic and alloy free.


Indications

ZirFit HT® Solid Zirconia is indicated for posterior crowns, bridges, and screw-retained or cement-retained implant crowns. It is an esthetic alternative to PFMs or full-cast restorations. ZirFit HT® restorations are also ideal for patients who don't have the adequate preparation space for a PFM.

Tech tips

To ensure the correct final shade of the restoration, please make sure to define the Dentin shade on your Rx form. Dark prepped areas can have a dramatic effect on the final shade of the restoration.

Posterior Prep



Anterior Prep

